

AZGP1 Antibody, FITC conjugated



PACO32946

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Stimulates lipid degradation in adipocytes and causes the extensive fat losses associated with some advanced cancers. May bind polyunsaturated fatty acid, .

Gene ID:

AZGP1

Uniprot

P25311

Synonyms:

Zinc-alpha-2-glycoprotein (Zn-alpha-2-GP) (Zn-alpha-2-glycoprotein), AZGP1, ZAG ZNGP1

Immunogen:

Recombinant Human Zinc-alpha-2-glycoprotein protein (61-210AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A