

KRT8 Antibody, Biotin conjugated



PACO32987

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Together with KRT19, helps to link the contractile apparatus to dystrophin at the costameres of striated muscle.

Gene ID:

KRT8

Uniprot

P05787

Synonyms:

Keratin, type II cytoskeletal 8 (Cytokeratin-8) (CK-8) (Keratin-8) (K8) (Type-II keratin Kb8), KRT8, CYK8

Immunogen:

Recombinant Human Keratin, type II cytoskeletal 8 protein (2-483AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A