

SERPINB3 Antibody, FITC conjugated



PACO32990

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May act as a papain-like cysteine protease inhibitor to modulate the host immune response against tumor cells. Also functions as an inhibitor of UV-induced apoptosis via suppression of the activity of c-Jun NH(2)-terminal kinase (JNK1).

Gene ID:

SERPINB3

Uniprot

P29508

Synonyms:

Serpin B3 (Protein T4-A) (Squamous cell carcinoma antigen 1) (SCCA-1), SERPINB3, SCCA SCCA1

Immunogen:

Recombinant Human Serpin B3 protein (2-383AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A