

# SERPINB3 Antibody, Biotin conjugated



PACO32991

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May act as a papain-like cysteine protease inhibitor to modulate the host immune response against tumor cells. Also functions as an inhibitor of UV-induced apoptosis via suppression of the activity of c-Jun NH(2)-terminal kinase (JNK1).

**Gene ID:**

SERPINB3

**Uniprot**

P29508

**Synonyms:**

Serpin B3 (Protein T4-A) (Squamous cell carcinoma antigen 1) (SCCA-1), SERPINB3, SCCA SCCA1

**Immunogen:**

Recombinant Human Serpin B3 protein (2-383AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A