

CSH1 Antibody, HRP conjugated



PACO33013

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Produced only during pregnancy and is involved in stimulating lactation, fetal growth and metabolism. Does not interact with GHR but only activates PRLR through zinc-induced dimerization.

Gene ID:

CSH1

Uniprot

P0DML2

Synonyms:

Chorionic somatomammotropin hormone 1 (Choriomammotropin) (Lactogen) (Placental lactogen) (PL), CSH1

Immunogen:

Recombinant Human Chorionic somatomammotropin hormone 1 protein (27-217AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A