

COL3A1 Antibody, Biotin conjugated



PACO33039

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	COL3A1
Source:	Uniprot
Rabbit	P02461
Isotype:	Synonyms:
IgG	Collagen alpha-1(III) chain, COL3A1
Applications:	Immunogen:
ELISA	Recombinant Human Collagen alpha-1(III) chain protein (1222-1465AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A