

NAPSA Antibody, FITC conjugated



PACO33062

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in processing of pneumocyte surfactant precursors.

Gene ID:

NAPSA

Uniprot

O96009

Synonyms:

Napsin-A (EC 3.4.23. -) (Aspartyl protease 4) (ASP4) (Asp 4) (Napsin-1) (TA01/TA02), NAPSA, NAP1 NAPA

Immunogen:

Recombinant Human Napsin-A protein (64-420AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A