

ABCB1 Antibody, HRP conjugated



PACO33105

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Energy-dependent efflux pump responsible for decreased drug accumulation in multidrug-resistant cells.

Gene ID:

ABCB1

Uniprot

P08183

Synonyms:

Multidrug resistance protein 1 (EC 3.6.3.44) (ATP-binding cassette sub-family B member 1) (P-glycoprotein 1) (CD antigen CD243), ABCB1, MDR1 PGY1

Immunogen:

Recombinant Human Multidrug resistance protein 1 protein (236-297AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A