

TSHB Antibody, Biotin conjugated



PACO33119

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Indispensable for the control of thyroid structure and metabolism.

Gene ID:

TSHB

Uniprot

P01222

Synonyms:

Thyrotropin subunit beta (Thyroid-stimulating hormone subunit beta) (TSH-B) (TSH-beta) (Thyrotropin beta chain) (Thyrotropin alfa), TSHB

Immunogen:

Recombinant Human Thyrotropin subunit β protein (21-132AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A