

MSTN Antibody, HRP conjugated



PACO33145

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Acts specifically as a negative regulator of skeletal muscle growth.

Gene ID:

MSTN

Uniprot

O14793

Synonyms:

Growth/differentiation factor 8 (GDF-8) (Myostatin), MSTN, GDF8

Immunogen:

Recombinant Human Growth/differentiation factor 8 protein (267-375AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A