

# MSTN Antibody, FITC conjugated



PACO33146

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Acts specifically as a negative regulator of skeletal muscle growth.

**Gene ID:**

MSTN

**Uniprot**

O14793

**Synonyms:**

Growth/differentiation factor 8 (GDF-8) (Myostatin), MSTN, GDF8

**Immunogen:**

Recombinant Human Growth/differentiation factor 8 protein (267-375AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A