

ALPI Antibody, FITC conjugated



PACO33158

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

ALPI

Uniprot

P09923

Synonyms:

Intestinal-type alkaline phosphatase (IAP) (Intestinal alkaline phosphatase) (EC 3.1.3.1), ALPI

Immunogen:

Recombinant Human Intestinal-type alkaline phosphatase protein (26-500AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A