

ALPL Antibody, FITC conjugated



PACO33162

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

This isozyme may play a role in skeletal mineralization.

Gene ID:

ALPL

Uniprot

P05186

Synonyms:

Alkaline phosphatase, tissue-nonspecific isozyme (AP-TNAP) (TNSALP) (EC 3.1.3.1) (Alkaline phosphatase liver/bone/kidney isozyme), ALPL

Immunogen:

Recombinant Human Alkaline phosphatase, tissue-nonspecific isozyme protein (28-501AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A