

ICA1 Antibody, HRP conjugated



PACO33173

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in neurotransmitter secretion.

Gene ID:

ICA1

Uniprot

Q05084

Synonyms:

Islet cell autoantigen 1 (69 kDa islet cell autoantigen) (ICA69) (Islet cell autoantigen p69) (ICAp69) (p69), ICA1

Immunogen:

Recombinant Human Islet cell autoantigen 1 protein (1-268AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A