

SPA17 Antibody, Biotin conjugated



PACO33179

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Sperm surface zona pellucida binding protein. Helps to bind spermatozoa to the zona pellucida with high affinity. Might function in binding zona pellucida and carbohydrates.

Gene ID:

SPA17

Uniprot

Q15506

Synonyms:

Sperm surface protein Sp17 (Cancer/testis antigen 22) (CT22) (Sp17-1) (Sperm autoantigenic protein 17) (Sperm protein 17), SPA17, SP17

Immunogen:

Recombinant Human Sperm surface protein Sp17 protein (1-146AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A