

# HBA1 Antibody



PACO33188

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Involved in oxygen transport from the lung to the various peripheral tissues.

**Gene ID:**

HBA1

**Uniprot**

P69905

**Synonyms:**

Hemoglobin subunit alpha-1 (Alpha-1-globin) (Hemoglobin alpha-1 chain), hba1

**Immunogen:**

Recombinant Human Hemoglobin subunit &alpha; protein (2-142AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A