

PYY Antibody, Biotin conjugated



PACO33199

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

PYY

Uniprot

P10082

Synonyms:

Peptide YY (PYY) (PYY-I) (Peptide tyrosine tyrosine) [Cleaved into: Peptide YY(3-36) (PYY-II)], PYY

Immunogen:

Recombinant Human Peptide YY protein (31-64AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A