

# TPO Antibody



PACO33204

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Iodination and coupling of the hormonogenic tyrosines in thyroglobulin to yield the thyroid hormones T3 and T4.

**Gene ID:**

TPO

**Uniprot**

P07202

**Synonyms:**

Thyroid peroxidase (TPO) (EC 1.11.1.8), TPO

**Immunogen:**


Recombinant Human Thyroid peroxidase protein (19-846AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human thyroid tissue using PACO33204 at dilution of 1:100.