

TPO Antibody, Biotin conjugated



PACO33207

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Iodination and coupling of the hormonogenic tyrosines in thyroglobulin to yield the thyroid hormones T3 and T4.

Gene ID:

TPO

Uniprot

P07202

Synonyms:

Thyroid peroxidase (TPO) (EC 1.11.1.8), TPO

Immunogen:

Recombinant Human Thyroid peroxidase protein (19-846AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A