

PRL Antibody, Biotin conjugated



PACO33227

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Prolactin acts primarily on the mammary gland by promoting lactation.

Gene ID:

PRL

Uniprot

P01236

Synonyms:

Prolactin (PRL), PRL

Immunogen:

Recombinant Human Prolactin protein (29-227AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A