

argF Antibody



PACO33280

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

argF

Uniprot

P06960

Synonyms:

Ornithine carbamoyltransferase subunit F (OTCase-2) (EC 2.1.3.3), argF

Immunogen:

Recombinant Escherichia coli Ornithine carbamoyltransferase subunit F protein (2-334AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A