

# CANT1 Antibody, Biotin conjugated



PACO33299

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CANT1

**Uniprot**

Q8WVQ1

**Synonyms:**

Soluble calcium-activated nucleotidase 1 (SCAN-1) (EC 3.6.1.6) (Apyrase homolog) (Putative MAPK-activating protein PM09) (Putative NF-kappa-B-activating protein 107), CANT1, SHAPY

**Immunogen:**

Recombinant Human Soluble calcium-activated nucleotidase 1 protein (80-397AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A