

MIB1 Antibody, Biotin conjugated



PACO33327

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

MIB1

Uniprot

Q86YT6

Synonyms:

E3 ubiquitin-protein ligase MIB1 (EC 2.3.2.27) (DAK-Interacting protein 1) (DIP-1) (Mind bomb homolog 1) (RING-type E3 ubiquitin transferase MIB1) (Zinc finger ZZ type with ankyrin repeat domain protein 2), MIB1, DIP1 KIAA1323 ZZANK2

Immunogen:

Recombinant Human E3 ubiquitin-protein ligase MIB1 protein (5-332AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A