

# ZMAT2 Antibody, FITC conjugated



PACO33354

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ZMAT2

**Uniprot**

Q96NC0

**Synonyms:**

Zinc finger matrin-type protein 2, ZMAT2

**Immunogen:**

Recombinant Human Zinc finger matrin-type protein 2 protein (1-194AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A