

TROVE2 Antibody, HRP conjugated



PACO33373

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

RNA-binding protein that binds to several small cytoplasmic RNA molecules known as Y RNAs. May stabilize these RNAs from degradation. May play roles in cilia formation and/or maintenance.

Gene ID:

TROVE2

Uniprot

P10155

Synonyms:

60 kDa SS-A/Ro ribonucleoprotein (60 kDa Ro protein) (60 kDa ribonucleoprotein Ro) (RoRNP) (Ro 60 kDa autoantigen) (Sjogren syndrome antigen A2) (Sjogren syndrome type A antigen) (SS-A) (TROVE domain family member 2), TROVE2, RO60 SSA2

Immunogen:

Recombinant Human 60 kDa SS-A/Ro ribonucleoprotein protein (3-535AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A