

ZSCAN21 Antibody



PACO33380

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Strong transcriptional activator.

Gene ID:

ZSCAN21

Uniprot

Q9Y5A6

Synonyms:

Zinc finger and SCAN domain-containing protein 21 (Renal carcinoma antigen NY-REN-21) (Zinc finger protein 38 homolog) (Zfp-38), ZSCAN21, ZFP38 ZNF38

Immunogen:

Recombinant Human Zinc finger and SCAN domain-containing protein 21 protein (1-461AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A