

# ZSCAN21 Antibody, Biotin conjugated



PACO33383

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Strong transcriptional activator.

**Gene ID:**

ZSCAN21

**Uniprot**

Q9Y5A6

**Synonyms:**

Zinc finger and SCAN domain-containing protein 21 (Renal carcinoma antigen NY-REN-21) (Zinc finger protein 38 homolog) (Zfp-38), ZSCAN21, ZFP38 ZNF38

**Immunogen:**

Recombinant Human Zinc finger and SCAN domain-containing protein 21 protein (1-461AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A