

# KRAS Antibody, Biotin conjugated



PACO33395

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Ras proteins bind GDP/GTP and possess intrinsic GTPase activity.

**Gene ID:**

KRAS

**Uniprot**

P01116

**Synonyms:**

GTPase KRas (K-Ras 2) (Ki-Ras) (c-K-ras) (c-Ki-ras) [Cleaved into: GTPase KRas, N-terminally processed], KRAS, KRAS2 RASK2

**Immunogen:**

Recombinant Human GTPase KRas protein (2-168AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A