

KLK10 Antibody



PACO33404

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Has a tumor-suppressor role for NES1 in breast and prostate cancer.

Gene ID:

KLK10

Uniprot

O43240

Synonyms:

Kallikrein-10 (EC 3.4.21. -) (Normal epithelial cell-specific 1) (Protease serine-like 1), KLK10, NES1 PRSSL1

Immunogen:

Recombinant Human Kallikrein-10 protein (31-274AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A