

AGER Antibody, FITC conjugated



PACO33418

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

AGER

Uniprot

Q15109

Synonyms:

Advanced glycosylation end product-specific receptor (Receptor for advanced glycosylation end products), AGER, RAGE

Immunogen:

Recombinant Human Advanced glycosylation end product-specific receptor protein (23-342AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A