

CYS4 Antibody



PACO33472

Product Information

Size:

50ug

Reactivity:

Saccharomyces cerevisiae

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

CYS4

Uniprot

P32582

Synonyms:

Cystathionine beta-synthase (EC 4.2.1.22) (Beta-thionase) (Serine sulfhydrase) (Sulfur transfer protein 4), CYS4, STR4

Immunogen:

Recombinant Saccharomyces cerevisiae Cystathionine beta-synthase protein (1-345AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A