

ORM2 Antibody, HRP conjugated



PACO33481

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Functions as transport protein in the blood stream. Binds various hydrophobic ligands in the interior of its beta-barrel domain. Also binds synthetic drugs and influences their distribution and availability. Appears to function in modulating the activity of the immune system during the acute-phase reaction.

Gene ID:

ORM2

Uniprot

P19652

Synonyms:

Alpha-1-acid, glycoprotein 2 (AGP 2) (Orosomuroid-2) (OMD 2), ORM2, AGP2

Immunogen:

Recombinant Human Alpha-1-acid, glycoprotein 2 protein (19-201AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A