

PACO33512

Product Information

Size:

50ug

Reactivity:

Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:**Gene ID:**

Egln3

Uniprot

Q91UZ4

Synonyms:

Egl nine homolog 3 (EC 1.14.11.29) (Hypoxia-inducible factor prolyl hydroxylase 3) (HIF-PH3) (HIF-prolyl hydroxylase 3) (HPH-3) (Prolyl hydroxylase domain-containing protein 3) (PHD3) (SM-20), Egln3


Immunogen:

Recombinant Mouse Egl nine homolog 3 protein (2-239AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western Blot

Positive WB detected in: Rat brain tissue, Mouse brain tissue

All lanes: Egln3 antibody at 3.6 μ g/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 28 kDa

Observed band size: 35 kDa

.