

HLA-DQA1 Antibody, Biotin conjugated



PACO33531

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

HLA-DQA1

Uniprot

P01909

Synonyms:

HLA class II histocompatibility antigen, DQ alpha 1 chain (DC-1 alpha chain) (DC-alpha) (HLA-DCA) (MHC class II DQA1), HLA-DQA1


Immunogen:

Recombinant Human HLA class II histocompatibility antigen, DQ α 1 chain protein (24-213AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO33531 at dilution of 1:100.