

Product Information

Size:

50ug

Reactivity:

Mycobacterium tuberculosis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

MPT63 is a small, major secreted protein of unknown function from Mycobacterium tuberculosis that has been shown to have immunogenic properties and has been implicated in virulence.

Gene ID:

mpt63

Uniprot

P9WIP1

Synonyms:

Immunogenic protein MPT63 (16 kDa immunoprotective extracellular protein) (Antigen MPT63), mpt63

Immunogen:

Recombinant Mycobacterium tuberculosis Immunogenic protein MPT63 protein (30-159AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A