

PACO33540

Product Information

Size:

50ug

Reactivity:

Mycobacterium tuberculosis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Part of the ABC transporter complex PstSACB involved in phosphate import.

Gene ID:

pstS1

Uniprot

P9WGU1

Synonyms:

Phosphate-binding protein PstS 1 (PBP 1) (PhoS1) (PstS-1) (38-kDa glycolipoprotein) (38-kDa lipoprotein) (P38) (Antigen Ag78) (Protein antigen B) (Pab), pstS1, phoS1

Immunogen:

Recombinant Mycobacterium tuberculosis Phosphate-binding protein PstS 1 protein (24-373AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A