

# mpt51 Antibody



PACO33544

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mycobacterium tuberculosis

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May have a role in host tissue attachment, whereby ligands may include the serum protein fibronectin and small sugars.

**Gene ID:**

mpt51

**Uniprot**

P9WQN7

**Synonyms:**

MPT51/MPB51 antigen, mpt51, fbpC1 fbpD mpb51

**Immunogen:**

Recombinant Mycobacterium tuberculosis MPT51/MPB51 antigen protein (27-299AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A