

CTSK Antibody, FITC conjugated



PACO33574

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Closely involved in osteoclastic bone resorption and may participate partially in the disorder of bone remodeling. Displays potent endoprotease activity against fibrinogen at acid, pH. May play an important role in extracellular matrix degradation.

Gene ID:

CTSK

Uniprot

P43235

Synonyms:

Cathepsin K (EC 3.4.22.38) (Cathepsin O) (Cathepsin O2) (Cathepsin X), CTSK, CTSO CTSO2

Immunogen:

Recombinant Human Cathepsin K protein (115-329AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A