

PACO33580

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Thiol protease which is believed to participate in intracellular degradation and turnover of proteins. Has also been implicated in tumor invasion and metastasis.

Gene ID:

CTSF

Uniprot

Q9UBX1

Synonyms:

Cathepsin F (CATSF) (EC 3.4.22.41), CTSF

Immunogen:

Recombinant Human Cathepsin F protein (273-484AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A