

Hras Antibody, Biotin conjugated



PACO33615

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Ras proteins bind GDP/GTP and possess intrinsic GTPase activity.

Gene ID:

Hras

Uniprot

P20171

Synonyms:

GTPase HRas (H-Ras-1) (Transforming protein p21) (c-H-ras) (p21ras) [Cleaved into: GTPase HRas, N-terminally processed], Hras, Hras1

Immunogen:

Recombinant Rat GTPase HRas protein (2-186AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A