

PACO33684

Product Information

Size:

50ug

Reactivity:

Baker's yeast

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Septum formation and repair, especially under certain adverse conditions.

Gene ID:

CHS1

Uniprot

P08004

Synonyms:

Chitin synthase 1 (EC 2.4.1.16) (Chitin-UDP acetyl-glucosaminyl transferase 1), CHS1

Immunogen:

Recombinant Baker's yeast Chitin synthase 1 protein (4-200AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A