

PACO33700

---

## Product Information

**Size:**

50ug

**Reactivity:**

Salmonella typhi

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Salmonella (S. ) typhi is an important intracellular pathogen. Among the more than 2,300 closely-related Salmonella serovars bacteria recognized, S.

**Gene ID:**

viaA

**Uniprot**

Q8Z2R0

**Synonyms:**

Protein ViaA (VWA domain protein interacting with AAA ATPase), viaA

**Immunogen:**

Recombinant Salmonella typhi Protein ViaA protein (3-476AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A