

SPI1 Antibody



PACO33712

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

SPI1

Uniprot

P17947

Synonyms:

Transcription factor PU.1 (31 kDa-transforming protein), SPI1

Immunogen:

Recombinant Human Transcription factor PU.1 protein (1-270AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A