

# clpS Antibody, HRP conjugated



PACO33729

---

## Product Information

**Size:**

50ug

**Reactivity:**

Escherichia coli

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Involved in the modulation of the specificity of the ClpAP-mediated ATP-dependent protein degradation.

**Gene ID:**

clpS

**Uniprot**

P0A8Q6

**Synonyms:**

ATP-dependent Clp protease adapter protein ClpS, clpS, yljA

**Immunogen:**

Recombinant Escherichia coli ATP-dependent Clp protease adapter protein ClpS protein (1-105AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A