

clpS Antibody, Biotin conjugated



PACO33731

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in the modulation of the specificity of the ClpAP-mediated ATP-dependent protein degradation.

Gene ID:

clpS

Uniprot

P0A8Q6

Synonyms:

ATP-dependent Clp protease adapter protein ClpS, clpS, yljA

Immunogen:

Recombinant Escherichia coli ATP-dependent Clp protease adapter protein ClpS protein (1-105AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A