

MT1E Antibody, HRP conjugated



PACO33801

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Metallothioneins have a high content of cysteine residues that bind various heavy metals; these proteins are transcriptionally regulated by both heavy metals and glucocorticoids.

Gene ID:

MT1E

Uniprot

P04732

Synonyms:

Metallothionein-1E (MT-1E) (Metallothionein-IE) (MT-IE), MT1E

Immunogen:

Recombinant Human Metallothionein-1E protein (4-59AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A