

MT4 Antibody



PACO33832

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Seems to bind zinc and copper. Could play a special role in regulating zinc metabolism during the differentiation of stratified epithelia.

Gene ID:

MT4

Uniprot

P47944

Synonyms:

Metallothionein-4 (MT-4) (Metallothionein-IV) (MT-IV), MT4

Immunogen:

Recombinant Human Metallothionein-4 protein (1-60AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A