

HLA-DRA Antibody, Biotin conjugated



PACO33871

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

HLA-DRA

Uniprot

P01903

Synonyms:

HLA class II histocompatibility antigen, DR alpha chain (MHC class II antigen DRA), HLA-DRA, HLA-DRA1

Immunogen:

Recombinant Human HLA class II histocompatibility antigen, DR α chain protein (28-254AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A