

# SLC39A1 Antibody, Biotin conjugated



PACO33879

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Mediates zinc uptake. May function as a major endogenous zinc uptake transporter in many cells of the body. Responsible for the rapid uptake and accumulation of physiologically effective zinc in prostate cells.

**Gene ID:**

SLC39A1

**Uniprot**

Q9NY26

**Synonyms:**

Zinc transporter ZIP1 (Solute carrier family 39 member 1) (Zinc-iron-regulated transporter-like) (Zrt- and Irt-like protein 1) (ZIP-1) (hZIP1), SLC39A1, IRT1 ZIP1 ZIRTL

**Immunogen:**

Recombinant Human Zinc transporter ZIP1 protein (126-179AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A