

NTH1 Antibody



PACO33908

Product Information

Size:

50ug

Reactivity:

Saccharomyces cerevisiae

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Purification and characterization of neutral trehalase from the yeast ABYS1 mutant.

Gene ID:

NTH1

Uniprot

P32356

Synonyms:

Neutral trehalase (EC 3.2.1.28) (Alpha, alpha-trehalase) (Alpha, alpha-trehalose glucohydrolase), NTH1, NTH

Immunogen:

Recombinant Saccharomyces cerevisiae Neutral trehalase protein (546-751AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A